

Operate With
Enlitic Intelligence.



Curie|ENDEX

DRIVE EFFICIENCY AND ENHANCE PRODUCTIVITY




Eliminate hanging protocol display issues
Reduce manual intervention for radiologists and PACS Administrators
Improve radiologist satisfaction with your PACS
Lessen potential repetitive stress injuries
Enrich data with reliable, clinically relevant labels and metadata
Enhance workflow orchestration



About
Enlitic

Founded in 2014, Enlitic is a medical intelligence company that believes in evolving patient care. Our products, innovations, insights and solutions were developed with Enlitic Intelligence - Intelligence designed for Collaboration, Adaptability, Relevance, Enlightenment. Or as we say, Intelligence that cares.

LEARN MORE AT [ENLITIC.COM](https://enlitic.com)



Curie|ENDEX

Medical Image Standardization

THE PROBLEM

In today's medical imaging department countless hours are wasted dealing with problems that affect hanging protocols, image routing, and AI orchestration. At the core of our healthcare interoperability challenges is non-standardized clinical data.

- Series are technically identical but labeled differently creating display and processing issues
- Series descriptions lack clinical relevance
- Different facilities label studies differently, making comparisons impossible
- Data is missing or inaccurate
- And more...

WHAT WE DO

Curie|ENDEX transforms medical imaging data to a consistent, clinically relevant standard nomenclature and enables relevant clinical content linkage across different IT systems. Endex enables consistent display of hanging protocols, improves image routing and orchestration intelligence.

4x

Reduced mouse
mileage

2-5x

Return on
investment

16s

Saved per
study

HOW IT WORKS

- Endex maps heterogeneous DICOM headers to a defined universal ontology
- Studies and series now appear with normalized descriptions, enabling hanging protocols to work
- Relevant, uniform descriptions facilitate faster reporting time, easily queried studies, and less frustration from radiologists
- Manual intervention is reduced by 15-90 seconds per study
- Admin staff reclaim time from mundane tasks

02230799133...	02230799133006...	Feb-08-2022 12:15 PM	CT CHEST W/IV CONTRAST	CT	36702265689920...
Basic Viewer					
Description	Series	Modality	Instances		
(empty)	2	CT	2		
Axial Body Std. Axial CE	4	CT	121		
Axial Lung Std. Axial CE	5	CT	361		
Coronal Body Std. Volume CE Coronal	7	CT	92		
Sagittal Body Std. Volume CE Sagittal	8	CT	120		
Ax-MIP CTA Body CE Ax-MIP	10	CT	181		
Cor-MIP CTA Body CE Cor-MIP	11	CT	58		
Sag-MIP CTA Body CE Sag-MIP	12	CT	73		

02230799133...	02230799133006...	Feb-08-2022 12:15 PM	CT CHEST SINGLE	CT	ST688892670104...
Basic Viewer					
Description	Series	Modality	Instances		
LOCALIZER	2	CT	2		
CHEST AXIAL C+ THICK	4	CT	121		
CHEST AXIAL C+ LUNG THIN	5	CT	361		
CHEST CORONAL C+ THICK	7	CT	92		
CHEST SAGITTAL C+ THICK	8	CT	120		
CHEST AXIAL C+ MIP	10	CT	181		
CHEST CORONAL C+ MIP	11	CT	58		
CHEST SAGITTAL C+ MIP	12	CT	73		

Curie|ENDEX WILL:

- **Enhance** your workflow and **increase** operational productivity
- **Improve** radiologist job satisfaction
- **Reduce** your physician burnout
- **Decrease** the complexities of PACS hanging protocols
- **Expand** capacity across multiple roles



3420 E Harmony Rd Suite 125 Fort Collins CO 80528
enlitic.com 877.365.4842

Book a Demo Today:

INFO@ENLITIC.COM

INTELLIGENCE THAT CARES